

Legionella MWY Selective Agar

Shorthand Symbol:	MWY
Item Number:	40-1604
Format:	Petri dish plates, 90mm
Colour:	Black
Storage:	Constant temperature, protected from light, at 4-25 °C
Shelf-life:	?
pH:	6,8 ± 0.2

Intended use and fields of application (Specification)

Isolation and cultivation of *Legionella Spec.*

Composition (g/l)

Activated carbon	2
Yeast extract	10
ACES Buffer	10
Potassium hydroxide	2.1
L-Cysteine	0.4
Iron pyrophosphate	0.25
α-ketoglutarate	1
Glycine	3
Vancomycin	0.001
Polymyxin B (IU)	50000
Anisomycin	0.08
Bromothymol Blue	0,01
Bromocresol Purple	0,01
Agar	13

Microbiological quality control

Xebios is certified according to DIN EN ISO 9001 as well as DIN EN ISO 13485. The microbiological quality control is carried out by our laboratory accredited according to DIN EN ISO / IEC 17025 (registration number of the certificate: D-PL-17997-01-00).